

Evaluation of a Small Field of View Personal Gamma Spectrometer Under ^{137}Cs Irradiation Conditions



K. Adamis, A. Metallinos, I. Kandarakis and **E. David**

*Radiation Physics, Materials Technology and Biomedical Imaging Laboratory
Department of Biomedical Engineering, University of West Attica, Athens, Greece*



**Sixth International Conference
ISMART 2018
Engineering of Scintillation Materials and
Radiation Technologies**

9 - 12 October 2018

Nuclear energy

2

- ✓ Raise of the usage of nuclear energy in various scientific sectors.



✓ More often nuclear accidents or leakages.

✓ Need for measurement of radioactivity in everyday life increased.

Personal Radiation Spectrometer Detectors

3

- A gadget that may interest more and more people over the years.
- Characteristics
 - ✓ Isotope identification
 - ✓ Compactness
 - ✓ Portability (size, weight, etc)
 - ✓ High Sensitivity
 - ✓ **Angular response**

Materials I

4

- ✓ ^{137}Cs radio source of $0.87\mu\text{Ci}$



- ✓ GAGG:Ce inorganic scintillator crystals



Furukawa Co., LTD	Epic Crystal
$3\times 3\times 4\text{ mm}^3$	$3\times 3\times 5\text{ mm}^3$
up to $3\times 3\times 20\text{ mm}^3$ (7 samples)	$3\times 3\times 10\text{ mm}^3$



- ✓ with all crystal surfaces polished.

Materials II

5

SiPM characteristics

Company name	KETEK
Model	PM3350 trench design
Pixel size	$3 \times 3 \text{ mm}^2$
Cell size	$50 \mu\text{m}$
Cells per pixel	3600
Bias Voltage	+30.5 V
Gain	10^6
Dark rate	$< 300 \text{ kHz/mm}^2$
Photon Detection Efficiency at 520nm	27 %



Materials III


6

- ✓ **DT5720 Caen digitizer**
 - ✓ Desktop waveform digitizer
 - ✓ 4 channels
 - ✓ 12-bit 250 MS/s
 - ✓ 2 V_{pp} single ended input dynamics
 - ✓ Adjustable DC offset via 16-bit DAC on each channel in the ± 1 V range
- ✓ **Digital pulse processing for charge integration (DPP-CI) Control Software**



Experimental Setup I

7





Experimental Setup II

8

- Angular detection response in 5 positions, 0, 45, 90, 135 & 180 degrees.
- 30' measurement at each angle position.


The 2 experimental setups for GAGG $3 \times 3 \times 5 \text{ mm}^3$ & GAGG:Ce $3 \times 3 \times 10 \text{ mm}^3$



Results and Discussion I

9

30.5 V bias was chosen, as the one that provides the lowest energy resolution values.



Results and Discussion II

10

Comparison between Furukawa & Epic Crystals irradiated with ^{137}Cs for 30 min. The source was placed just upon the crystal for the excitation.

Crystal Thickness (mm)	Furukawa Crystals	Epic Crystals
	R (%)	R (%)
4	7.2	-
5	6.8	5.8
6	8.9	-
8	9.4	-
10	10.4	9.5
15	6	-
20	10.2	-


Sensitivity

11

Higher sensitivity for the bigger crystal

$$\text{Sensitivity} = \frac{\text{total counts}}{\text{emitted photons}}$$



Sensitivity response of Epic crystals



Energy Resolution (ER)

12

Energy resolution response of Epic crystals



MEAN VALUES

@5mm = 6.7%

@10mm = 9.4%.

Photo percentage

13

$$PP = \frac{662\text{keV photopeak counts}}{\text{total recorded counts}} \times 100\%$$

MEAN VALUES for ^{137}Cs

- 11.2% -5mm thick
- 15.1% -10mm thick



Angles (degrees)	Photofraction (%)	
	3x3x5mm ³	3x3x10mm ³
0	10.3	14.4
45	12.8	15.4
90	10.6	14.8
135	11.3	15.1
157	11.8	16.1

MEAN VALUES for ^{22}Na

- 13% -5mm thick
- 17% -10mm thick

Stability-in-time & Dark counts measurements

14



Future work....

15

- Additional features, such as
 1. detector response to other gamma energies,
 2. possible ability to construct a new DAC in order to communicate via WiFi with mobile devices via web browser,
 3. construct smaller electronics to have the ability to take place in a miniature device.

can also be implemented to improve this device.

Acknowledgement

16

«This research is implemented through IKY scholarships programme and co-financed by the European Union (European Social Fund - ESF) and Greek national funds through the action entitled "Reinforcement of Postdoctoral Researchers", in the framework of the Operational Programme "Human Resources Development Program, Education and Lifelong Learning" of the National Strategic Reference Framework (NSRF) 2014 – 2020».

References

- David S, Michail C, Valais I, Kalyvas N, Bakas A, Gektin A, Kandarakis I and Kourkoutas K 2015 *Proceedings of the SCinTE 2015 international conference, Athens, Greece, November 5-7*; ISBN: 978-960-98739-8-7
- David S, Georgiou M, Fysikopoulos E, and Loudos G 2015 *Phys. Med.* **31** 763-766
- Lipovec A, Shimazoe K and Takahashi H, 2016 *17th International Workshop on Radiation Imaging Detectors, JINST* **11** C03026 doi:10.1088/1748-0221/11/03/C03026
- Seitz B, Campos Rivera N and Stewart A G 2016 *IEEE Trans Nucl. Sci.* **63** 503-504
- KETEK GmbH, PM3350, 3x3 mm² - 50 μm microcell type, datasheet, available online at: <https://www.ketek.net/wp-content/uploads/2017/01/KETEK-PM3325-EB-PM3350-EB-Datasheet.pdf>.
- Stratos D, Georgiou M, Fysikopoulos E and Loudos G 2013 *Nucl. Instr. and Meth. Phys. Res. A* **702** 121–125
- Park H M and Joo K S 2016 *Sensor* **16** 919 DOI:[10.3390/s16060919](https://doi.org/10.3390/s16060919)
- Caen SpA, “MOD. DT5720, 4 channel 12 bit 250 MS/s Digitizer, Technical Information Manual”, Revision no.8, p.6-13, 2013
- Caen SpA, “UM2089, DPCI User Manual”, Rev.0-08, p.5-9, 2013
- Delimani G, Georgiou M, Papadimitroulas P, Papamichalis P, Mikropoulos K, David S and Loudos G 2014 *EJMP* **30** Supplement e110–e111 <http://dx.doi.org/10.1016/j.ejmp.2014.07.314>

Thank you for your attention